

INSTECH Flash Table-Top Autoclave B-Type

INSTECH Flash Table-Top Autoclave B-Type conform to the latest world standards.



Model: ISFB-25



Model: ISFB-25

- ### Product Description
- ✓ INSTECH Table-Top Autoclave ideal for Sterile applications in Hospitals, Medical Colleges, Research Centers & Pharma Units.
 - ✓ Our manufacturing process entails best-in-class engineering practices like CNC Laser Cutting & bending followed by TIG & Mechanised Finishing
 - ✓ Our products undergo a comprehensive QA regime and come with all safety devices required for a trusted and safe operation
 - ✓ We conduct thorough testing of parts for durability

- ### Salient Features
- ✓ B-Process, Vacuum Drying Process,
 - ✓ Audio Visual Alarm Notifies Temp., Pressure Deviation Pressure Sensor ,Power & Door Improper Closed Failures.
 - ✓ Built-in Open type Water Tank.
 - ✓ Vacuum drying System Diaphragm VACUUM PUMP.
 - ✓ Automatic Control With LCD Display.
 - ✓ Relief Valve in Case of High Pressure.
 - ✓ Performs B&D test, VLT & Helix test

Bespoke Customized Capabilities

- High Efficiency Ultimate Vacuum
- Sterilization Temp. (121°C- 134°C)
- Standard USB Ports/HMI
- Working Pressure(-0.95 - 2.5 bar)
- SS Tray Shelves 3 Nos.
- These all are available optionally

MODELS

Model	ISFB-25
Cap.	25 Ltr.
Dimensions (mm)	270Dia X 450Depth
Power Consumption	1950W