

INSTECH Ultrasonic Cleaner conform to the latest world standards.



Capacity 30 Ltr.



Capacity 20 Ltr.

Product Description

- ✓ INSTECH ULTRASONIC CLEANER is ideal for Cleaning applications in Hospitals, Medical Colleges, Research Centers & Pharma Units
- ✓ Use of Best in class manufacturing processes & engineering practices like CNC Laser Cutting & bending, TIG Welding & Mechanised Finishing
- ✓ Our products undergo a comprehensive QA regime and come with all safety devices required for a trusted and safe operation
- ✓ We conduct thorough testing of parts, thus ensuring reliability and durability

Salient Features

- ✓ Adjustable Digital Timer and Temperature Indicator
- ✓ Highly efficient cleaning can be done using industrial alcohol at High frequency
- ✓ A well built & reliable controller circuit
- ✓ Can be used in the Biochemistry, Microbiology
- ✓ Cleaning to be done using RO water only

Models

Models Name	Capacity
IS USC-20	20 Ltr.
IS USC-30	30 Ltr.
IS USC-40	40 Ltr.
IS USC-50	50 Ltr.
IS USC-80	80 Ltr.

Technical Features

Ultrasonic Freq.	40 KHz
Input Supply	AC 220-240V
Freq. in Hz	50 Hz.
Time Range	1-15 min
Ultrasonic Power	500W
Temp Range	30°C to 65°C
Heating Power	500/1000W